

ProVetScan AS

High quality of control in porcine,
ovine and caprine



Provetscan AS Wi-Fi Probe

High Quality Ultrasound scanner for reproductive diagnostic in Porcine, Ovine and Caprine.



Specifications:

Monitor: Smartphones and tablets of Android and iOS.

Frequency: 3,5 - 5 MHz

Scanning mode: Convex

Number of crystals: 80 crystals

Scanning depth: 90-180 mm




Software: Gain, Depth, Enhancement, Dynamic Range, Harmonic, save image and videos.

Grey levels: 256 levels

Battery working time: 5 hours

Size: 160x50x25 mm

Weight: 300 g



Included:

Provetscan AS Wi-Fi probe
Tablet Samsung Galaxy Tab A 8"
Protection cover
Neck belt
Transport suitcase

Optional:

Wireless charger

New VeTec S.L.