

ProVetScan SF-2

Back Fat and Loin diameter
measurement on P2



Provetscan SF-2 Wi-Fi Probe



High quality and accuracy for back fat and loin diameter in porcine and bovine

Specifications:

Monitor: Tablets of Android and iOS.

Frequency: 5 MHz

Scanning mode: Linear

Number of crystals: 80 crystals





Scanning depth: 45-120 mm

Software: Gain, Depth, back fat and loin measurement, save images and data.

Grey scales: 256 levels Battery working time: 3 horas

Main unit size: 155x60x20 mm

Weight: 300 g



Incluye:

Provetscan SF-2 Wi-Fi probe
Galaxy Tab A 8" or iPad mini
Protection cover
Neck belt
Transport suitcase

Optional:

Wireless charger

New VeTec S.L.