

# ProVetScan CR

Rectal convex reproductive diagnosis in  
cattle and horses



# Provetscan CR Wi-Fi Probe



Ultrasound scanner for rectal convex reproductive diagnostics in cows and mares.

## Specifications:

Monitor: iPad series.  
Frequency: 3,5-5 MHz  
Scanning mode: Convex  
Number of crystals: 80  
Scanning depth: 100-190 mm  
Display mode: Modo B  
Chip Wi-Fi: 2,4-5G  
Software: Gain, Depth, Enhancement, Dynamic Range, Cine-loop, Grid, measurements, gestation age calculation, save images and videos.  
Grey levels: 256  
Battery working time: 6 hours  
Main unit size: 180x55x23 mm  
Cable length: 1,80 m  
Weight: 400 g



VetUS



## Included:

Provetscan CR Wi-Fi probe  
iPad mini 64Gb  
Waterproof cover IP68  
Neck harness  
2A Charger  
Transpor Suitcase

## Optional:

Wireless charger

**New Ve Tec S.L.**