

# ProVetScan AS

High quality of control in porcine,  
ovine and caprine



# Provetscan AS Wi-Fi Probe



High Quality Ultrasound scanner for reproductive diagnostic in Porcine, Ovine and Caprine.

## Specifications:

Monitor: Smartphones and tablets of Android and iOS.

Frequency: 3,5 - 5 MHz

Scanning mode: Convex

Number of crystals: 80 crystals

Scanning depth: 90-180 mm

Software: Gain, Depth, Enhancement,




Dynamic Range, Harmonic, save image and videos.

Grey levels: 256 levels

Battery working time: 5 hours

Size: 160x50x25 mm

Weight: 300 g



## Included:

- Provetscan AS Wi-Fi probe
- Neck belt
- 2A Charger
- Wireless charger
- Transport suitcase

## Optional:

- iPad mini or Android Tablet
- Protection cover

**New VeTec S.L.**