

ProVetScan SF-2

Back Fat and Loin diameter
measurement on P2







Provetscan SF-2 Wi-Fi Probe



High quality and accuracy for back fat and loin diameter in porcine and bovine

Specifications:

Monitor: Tablets of Android and iOS.
Frequency: 5 MHz
Scanning mode: Linear
Number of crystals: 80 crystals
Scanning depth: 45-120 mm
Software: Gain, Depth, back fat and loin measurement, save images and data.
Grey scales: 256 levels
Battery working time: 3 hours
Main unit size: 155x60x20 mm
Weight: 300 g



Included:

Provetscan SF-2 Wi-Fi Probe
Neck belt
2A Charger
Wireless charger
Transport suitcase

Optional:

iPad mini or Android Tablet
Protection cover

New VeTec S.L.