

# ProVetScan SL-2S

High frequency and resolution  
wifi ultrasound probe



# Provetscan SL-2S Wi-Fi Probe



Best image quality for high frequency short distance diagnostics

## Specifications:

Monitor: Smartphones and tablets of Android and iOS.

Frequency: 7,5 - 10 MHz

Scanning mode: Linear

Number of crystals: 128 crystals

Scanning depth: 20-55 mm


Software: Gain, Depth, Enhancement, Dynamic Range, Harmonic, save images and videos.

Grey scales: 256 levels

Battery working time: 3 hours

Main unit size: 155x67x20 mm

Weight: 300 g



## Included:

Provetscan SL-2S Wi-Fi Probe

Neck belt

2A Charger

Wireless charger

Transport suitcase

## Optional:

iPad mini or Android Tablet

Protection cover



**New VeTec S.L.**