

ProVetScan SR-2C

High resolution ultrasound scanner for
diagnosis in cattle and horses



HD GLASSES





Provetscan SR-2C Wi-Fi Probe



Ultrasound scanner for reproductive diagnostics in cows, mares and pneumonia in calves.

Specifications:

Monitor: iPad series.
Frequency: 6,5-8 MHz
Scanning mode: Linear
Number of crystals: 128
Scanning depth: 40-100 mm
Display mode: Modo B y BM
Doppler color: Color, PDI, PW
Chip Wi-Fi: 2,4-5G
Software: Gain, Depth, Enhancement, Dynamic Range, Cine-loop, Grid, measurements, gestation age calculation, save images and videos.
Grey levels: 256
Battery working time: 6 hours (14 hours at Stand By)
Main unit size: 180x55x23 mm
Cable length: 1,80 m
Weight: 400 g



VetUS



Now with high-resolution HD USB-C Goggles to watch the best image definition:



Included:

Provetscan SR-2C Wi-Fi probe
Neck harness
2A Charger
Wireless charger
Transport suitcase

Optional:

iPad mini
Waterproof cover IP68
HD USB-C Goggles

New VeTec S.L.