

ProVetScan SS-3SA

The best image quality for
abdominal diagnosis



Provetscan SS-3SA Wi-Fi Probe



The best image quality for abdominal diagnosis in cows

Specifications:

Monitor: Smartphones and tablets of Android and iOS.

Frequency: 3,5 - 5 MHz

Scanning mode: Cónvex

Number of crystals: 128 crystals

Scanning depth: 90-305 mm

Software: Gain, Depth, Enhacement,


Dynamic Range, Harmonic, save images and videos.

Grey scales: 256 levels


Battery working time: 3 hours

Main unit size: 155x67x20 mm

Weight: 300 g



VetUS



Included:

- Provetscan SS-3SA Wi-Fi Probe
- Neck belt
- 2A Charger
- Wireless charger
- Transport suitcase

Optional:

- iPad mini or Android Tablet
- Protection cover



New VeTec S.L.