

# ProVetScan SU

Back Fat and Loin diameter and area measurement







# Provetscan SU-3A Wi-Fi Probe



High quality and accuracy for back fat and loin diameter and area in porcine and bovine

## Specifications:

Monitor: Android e iOS tablets.  
Frequency: 3,5 MHz  
Scanning mode: Linear  
Number of crystals: 256 crystals  
Scanning depth: 70-200 mm  
Software: Gain, Depth, back fat and loin, save images and data.  
Grey scales: 256 levels  
Battery working time: 3 hours  
Main unit size: 180x85x25 mm  
Weight: 400 g



## Include:

Provetscan SU-3A Wi-Fi probe  
Neck belt  
2A Charger  
Transport suitcase

## Optional:

Android Tablet or iPad mini  
Protection cover

**New VeTec S.L.**