

# ProVetScan SV-3

Reproductive control in porcine,  
ovine and caprine




# Provetscan SV-3 Wi-Fi Probe



Ultrasound scanner for reproductive diagnostic in Porcine, Ovine and Caprine.

## Specifications:

Monitor: Smartphones and tablets of Android and iOS.  
Frequency: 3,5 MHz  
Scanning mode: Sectorial  
Number of crystals: One annular piezo crystal  
Scanning depth: 100-180 mm  
Wi-Fi chip: 2,4 - 5G  
Grey levels: 256 niveles  
Battery working time: 5 hours  
Size: 140x45x40 mm  
Weight: 200 g



## Included:

Provetscan SV-3 Wi-Fi probe  
Neck belt  
2A Charger  
Wireless charger  
Transport suitcase

## Optional:

iPad mini or Android  
Tablet Protection cover

**New VeTec S.L.**